

OUR HEART AND SOUL...



FUMALEN[®]

SUTURES
FOR HEART VALVES

FUMEDICA 

A perfect tool results in a perfect outcome

Innovative Sutures from FUMEDICA

Advantages

- Easy removal due to optimized packaging
- Premium needle with enhanced surface
- Oval teflon pledgets
- Ideal gliding properties and knot stability due to polyester re-engineering
- Either coated (**Coated**)
- Or non coated (**Non coated**)



Basic product line

Needle bending	Needle length mm	USP	Suture length cm	Pledget / Needletype	Packaging	Order number	Order number
H	17	2-0	75	2	3	H17-2723 C	H17-2723 N
H	17	2-0	75	2	4	H17-2724 C	H17-2724 N
H	17	2-0	75	2	8	H17-2728 C	H17-2728 N
H	26	2-0	90	2	3	H26-0923 C	H26-0923 N
H	26	2-0	90	2	4	H26-0924 C	H26-0924 N
H	26	2-0	90	2	8	H26-0928 C	H26-0928 N

The shown sizes and combinations are an extract of the available product range. Upon request we are able to produce every needle and suture combination according to individual requirements.

Order number

Needle bending	Needle length mm	USP	Suture length cm	Pledget / Needletype	Packaging	Sutures type
H = ½ circle HD = ¾ circle	17 = 17mm 20 = 20mm 26 = 26mm	0 = 0 2 = 2-0 3 = 3-0	7 = 75 cm 9 = 90 cm	0 = taper point w/o pledget 1 = taper cut w/o pledget 2 = taper point with pledget 3 x 6mm 3 = taper cut with pledget 3 x 6mm	3 = 36 dozen single sutures green 4 = 10 pouches white 4 sutures 2 x green, 2 x white 8 = 10 pouches white 8 sutures 4 x green, 4 x white	Coated Non coated

Germany / International

Fumedica Medizintechnik GmbH
Ziegelwasen 4
D-72336 Balingen
Phone +49 (0)7471 920 00 50
Fax +49 (0)7471 920 00 52
Internet www.fumedica.de
E-Mail fumedica@fumedica.de

Switzerland

Fumedica AG
Luzernerstrasse 91
CH-5630 Muri
Phone +41 (0)56 675 91 00
Fax +41 (0)56 675 91 09
Internet www.fumedica.ch
E-Mail fumedica@fumedica.ch

