


OUR HEART AND SOUL...



STERNOTOMY
DOUBLE WIRES
MADE OF
STAINLESS STEEL
AND TITANIUM

A perfect tool results in a perfect outcome

Innovative Sternotomy Wires from FUMEDICA



Advantages

- High degree of break resistance at twisting point
- Large contact area at sternum
- Ease of use due to twister

Double Wires made of stainless steel, needle DS 48, sterile packed

Wire size USP	Wire length cm	Wires per pouch	No. pouches per box	Order number
7	2 x 30	1	12	DS4-770-1
7	2 x 30	2	12	DS4-770-2
7	2 x 30	3	12	DS4-770-3

Double Wires made of titanium, needle DS 48, sterile packed

Wire size USP	Wire length cm	Wires per pouch	No. pouches per box	Order number
7	2 x 30	1	12	TDS4-770-1
7	2 x 30	2	12	TDS4-770-2
7	2 x 30	3	12	TDS4-770-3

Twister device for twisting the double wire



Order number
S4-100

Germany / International

Fumedica Medizintechnik GmbH
Ziegelwasen 4
D-72336 Balingen
Phone +49 (0)7471 920 00 50
Fax +49 (0)7471 920 00 52
Internet www.fumedica.de
E-Mail fumedica@fumedica.de

Switzerland

Fumedica AG
Luzernerstrasse 91
CH-5630 Muri
Phone +41 (0)56 675 91 00
Fax +41 (0)56 675 91 09
Internet www.fumedica.ch
E-Mail fumedica@fumedica.ch

