


OUR HEART AND SOUL...




PREMIUM
STERNOTOMY
WIRES

A perfect tool results in a perfect outcome

Innovative Sternotomy Wires from FUMEDICA

The needle...

... has a special borehole and thus an optimized connection and a better border flexibility



The «Strong Wire»...

... has a significantly larger surface area and tensile strength and minimizes the risk of cutting into the bone structure

The «flexible Wire»...

...is 30 % more flexible than standard wire*

- Enhanced bending property
- Distance-free cerclage
- less sternum dehiscence

* compared to a normal wire with the same diameter



Special Wires extra flexible, needle FBX 48, sterile packed

Wire size USP	Wire lenght* cm	Wires per pouch	No. pouches per box	Order number
7	45	1	12	KS1-745-1
7	45	4	12	KS1-745-4
7	45	5	12	KS1-745-5
7	45	6	12	KS1-745-6
7	45	8	12	KS1-745-8

Special Wires extra flexible, needle FB 48, sterile packed

9	45	1	12	KDS4-945-1
9	45	4	12	KDS4-945-4

*Upon request we are able to produce other wire lengths, needle types or packaging variations.

Germany / International

Fumedica Medizintechnik GmbH
Ziegelwasen 4
D-72336 Balingen
Phone +49 (0)7471 920 00 50
Fax +49 (0)7471 920 00 52
Internet www.fumedica.de
E-Mail fumedica@fumedica.de

Switzerland

Fumedica AG
Luzernerstrasse 91
CH-5630 Muri
Phone +41 (0)56 675 91 00
Fax +41 (0)56 675 91 09
Internet www.fumedica.ch
E-Mail fumedica@fumedica.ch

